
















**Specimen Guide for Collections
GI Specimens**

Test Ordered	Inpatient/ED	Outpatient/Clinic/Urgent Care
C. difficile PCR	<p>Stool in screw cap container</p> 	<p>Stool in screw cap container</p> 
GI Panel	<p>Stool in screw cap container</p> 	<p>Fecal Swab</p>  <p>& Stool in screw cap container</p> 
<p>VRE Culture ESBL Culture Rectal Group A Strep Culture</p>	<p>ESwab</p> 	<p>ESwab</p> 

<p>Fecal Lactoferrin Calprotectin Fecal Fat Pancreatic Elastase</p>	<p>Stool in screw cap container</p> 	<p>Stool in screw cap container</p> 
<p>Fecal Ova and Parasite (OAP)</p>	<p>Stool in screw cap container</p> 	<p>EcoFix O&P</p> 
<p>H Pylori</p>	<p>Stool in screw cap container</p> 	<p>C & S Media</p> 
<p>Fecal Occult Blood</p>	<p>Stool in screw cap container</p> 	<p>OC-Light Sampling Bottle</p> 

Please see [Laboratory Specimen Labeling](#) policy for instructions on labeling specimens.